



No. Store/Purchase/Virology2/2015-AIIMS- RPR
2015

Dated: 25th February,

Inviting quotations for Virology for Microbiology, AIIMS, Raipur.

QUOTATION NOTICE

Sealed quotations are invited from intending manufacturer/stockist/distributors/ agents having registered TIN and relevant documents for Purchase of items for molecular diagnostics for **Virology** for department of **Microbiology**, AIIMS Raipur up to 3:00 pm of 03rd March, 2015. The quotations will be opened on the same day as per detailed below:-

S. No.	Items name	Specification	Probable vendor	Qty	Unit Rate in Rs.)
1.	Gloves, (7No.)	Sterile (Nitrile)	Tarson/	2 Pack x 100no.	
2.	Powder free sterile gloves 7No.	sterile disposable surgical rubber glove	Eppendorf	2 Pack x 100no.	
3.	SS Box(1Litre, 500ml)	Stainless steel	Tarson/	2 Box each	
4.	Rack to keep tubes of 1.5ml, 15ml tubes and 200ul capacity	Made of Polypropylene	Tarson/ Eppendorf	05 each	
5.	Conical centrifuge tubes with lid of 15ml capacity	Leak proof Polystyrene tubes with printed volume graduation and writing area along the length of the tube. Certified DNase, RNase & Pyrogen free	Tarson/ Eppendorf/ P ² fact	2 Pack x 500no.	
6.	Microcentrifuge tube with lid 1.5ml 200ul	Polypropylene, autoclavable	Tarson / Eppendorf/ Invitrogen/ Life technologies	2 pack x 1000No. 1 Pack of 1000 No.	
7.	Vacuum pump with tubing	Compact light aluminium die construction, oil free, piston powered, moisture trap/vacuum regulator, use of liquid filtration for microbiology.	Tarson/ Eppendorf/ Invitrogen/ Life technologies	1No.	
8.	Microtip 10ul 200ul 1000ul	Certified DNase, RNase & Pyrogen free	Eppendorf/ Tarson/ BR Biochem	10 pack each x 1000	
9.	Parafilm (1 roll of 4" X 125')	Thermoplastic self sealing film that holds moisture loss to a minimum and offers excellent barrier protection.	Tarson/ Abdos	1 roll of 4" x 125'	

10.	RNA extraction kit	Isolate RNA from body fluids	Qiagen	2x50 test	
11.	Hanks Balanced salt solution with phenol red indicator (10X)	Required for VTM preparation	Life technologies	1X 1000ml	
12.	Penicillin & Streptomycin Penicillin 10000U & 10mg Streptomycin/ml	Required for VTM preparation	/ Hi media	5X 50 ml	
13.	NaHCO ₃ (Sodium Bicarbonate)	Required for VTM preparation		1 X 100gm	
14.	Vortex	Continuous & Touch mode. Heavy metal base & rubber feet. Maximum speed 3000rpm	Tarson/ Eppendorf	01	

1. Taxes, if any (Kindly Mention).
2. Delivery Schedule – within 10 days from the date of issue of PO.
3. Price should be for Destination basis.
4. LD @ 0.5% of PO value per week for delay of supply of material subject to maximum up to 10% of PO value.

(Store Officer)
AIIMS, Raipur (CG)