



अखिल भारतीय आयुर्विज्ञान संस्थान, रायपुर(छत्तीसगढ़)
All India Institute of Medical Sciences, Raipur (Chhattisgarh)
Tatibandh, GE Road,
Raipur-492 099 (CG)
www.aiimsraipur.edu.in

Admin/ Store/Purchase/Non-Consumable/Pathology & Lab Medicine/2016-AIIMS.RPR
Date: 27.04.2016

Inviting Quotations for Purchase of Tube Roller & Micropipettes non-Consumable items for department of Pathology & Lab Medicine, AIIMS Raipur.


QUOTATION NOTICE

Sealed quotations are invited from intending registered Stockist / Distributors having TIN and relevant documents for Purchase of Tube Roller & Micropipettes non-Consumable items for department of Pathology & Lab Medicine, AIIMS Raipur. The quotation with copy of certificate of TIN & other documents should be submitted to office of **Store Officer** up to 05/05/2016 before 3:00 pm. The quotations will be opened on the same day at 3:30pm. Details of item are given as under:-

S.No.	Complete Description of items	Quantity Required	Purpose	UNIT RATE IN Rs.
1	TUBE ROLLER- 10 rollers <ul style="list-style-type: none">• Should have brushless DC motor for maintenance free long life.• Should have microprocessor controlled unit.• Should have 10 rollers.• Should have timer setting from 1 minute to 99 minutes & continuous mode• Should have easily removable rollers.• Should have Gentle Rocking & Rolling motion for optimum mixing.• Should have pulse mode function-momentary pause for optimum mixing.• Should have adjustable speed ranging from 20 to 60 RPM or better.• Should have max loading Capacity of 2Kg- for different application.• Should Comply with IP 21 protection class.• Should be able to use in incubators up to 40 °C and in cold rooms down to 5°C.• Company should provide 2 years warranty.	01 Set	For use in Haematology Lab in the Department of Pathology & Lab Medicine.	

2	<p>MICROPIPETTES:</p> <p>*Micropipette. Adjustable volume, fully autoclavable, Spring Loaded Tip Cone for connecting tips very tightly.</p> <p>*Adjustable opening for adjusting pipettes to a specific liquid and volume.</p> <p>*Control Button with very low operating force. Color indication for pipette volume.</p> <p>*Tip ejector with very low operating force, positioned for perfect ergonomics.</p> <p>*Volume Display : 4 Digits.</p> <p>*Piston should be manufactured in such material to provide thermal, mechanical and chemical stability.</p> <p>*Very easy removable lower part for cleaning pipette.</p> <p>*No discoloration upon UV irradiation</p> <p>*3 Years Warranty.</p> <p>*Company should provide 3 years calibration & repairing.</p> <p>*Volume range- 0.1-2.5 µl, 2.0 – 20 µl, 5 – 50 µl, 20 – 200 µl micro litres.</p> <p>*Should be supplied with individual quality certificate (QC) and calibration report according to ISO8655.</p> <ul style="list-style-type: none"> • Should have CE and IVD compliant. • Pipette stand, with 6 Pipettes should be provided. • Microtips- 10 µl (can take up to 20 µl) low retention with rack – 10 nos. x 96 tips should be provided. • Microtips- 200 µl low retention with rack – 10 nos x 96 tips should be provided. • Microtips- 1000 µl with rack – 06 nos x 96 tips should be provided. 	02 Kit/Set	For use in histopathology . Cytopathology. flow cytometry & hematology lab in Department of Pathology & Lab Medicine.
---	--	------------	---

- 1 Taxes, if any (Kindly mention the rate.)
- 2 Delivery Schedule - Within 15 to 21 days from the date of issue of PO
- 3 Price should be for- Destination basis
- 4 L.D @ 0.5% of delayed supply per week or part of week for delay of supply of material subject to maximum up to 10% of delayed supply should be deducted.
- 5 Quotation Name/No. should be mentioned on top of envelopes.



Sr. Administrative Officer
AIIMS, Raipur (C.G.)